ENGLISH WARNINGS

The Foam Roller is intended for commercial or home use. We strive to make the Foam Roller as safe as possible. This is an advanced mechanical device with electric components. If the Foam Roller and its accessories are not used or maintained properly, there is a risk of fire, electric shock, or injury. Do not use the Foam Roller on your head or near your genitals. Do not use the Foam Roller above the waist during pregnancy. Do not get hair caught in any moving part of the Foam Roller. Do not use the Foam Roller above your Adam's apple or C4. Do not get the Foam Roller wet. Only clean the Foam Roller with a magic eraser or a lightly damp towel. Do not block the vents of the motor. Never attempt to disassemble the Foam Roller. Never store the Foam Roller at temperatures above 40°C/104°F. Only charge the Foam Roller when the ambient temperature is between 0°C/32°F and 40°C/104°F. Never expose the Foam Roller to water. For long-term storage, the Foam Roller should be fully charged. This product is not intended to diagnose, treat, or prevent any diseases. BEFORE USING OR CHARGING THE Foam Roller, READ ALL THE INSTRUCTIONS AND CAUTIONARY MARKINGS IN THIS MANUAL AND ON THE Foam Roller. When using the Foam Roller, these basic precautions should always be followed including:

1. USE ONLY AS INSTRUCTED. Use the Foam Roller as described in this Foam Roller Instruction Manual. Do not carry out any maintenance other than what is shown in the Foam Roller Instruction Manual or advised by Flexir Recovery.
2. NOT FOR CHILDREN. The Foam Roller is not intended for use by young children or persons with reduced physical, sensory or reasoning capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction by a responsible person. Do not allow the Foam Roller to be used as a toy. Close attention is necessary when used by

or near children. Children should be supervised to ensure that they do not play with the Foam Roller.

3. **REDUCING RISK OF INJURY.** Charge the Foam Roller by plugging cable included in the package on your USB AC adapter of choice. Foam Roller can also be charged with a regular USB-C charger certified by a competent body. Do not use an AC adapter that is not certified or that has received a sharp blow, been dropped, or damaged in any way.

4. **CHARGING LOCATIONS.** The Foam Roller should be charged indoors in a well ventilated, dry location. Do not charge the Foam Roller outdoors, in a bathroom or within 10 feet (3.1 meters) of a bathtub, or pool. Do not use the Foam Roller on wet surfaces and do not expose the device to moisture, rain, or snow. Do not use the Foam Roller in the presence of explosive atmospheres (gaseous fumes, dust, or flammable materials). Sparks may be generated, possibly causing a fire.

5. **DO NOT OVERCHARGE.** Do not leave the Foam Roller connected to the adapter for more than an hour after the Foam Roller has been fully charged. The Foam Roller includes a protection system to avoid the risk of overcharging. However, overcharging the Foam Roller may reduce its life over time.

6. **DO NOT BURN OR INCINERATE THE Foam Roller.** The internal battery may explode, causing personal injury or damage. Toxic fumes and materials are created when the battery is burned.

7. **DO NOT CRUSH, DROP, OR DAMAGE THE Foam Roller charging cable.** Do not use a charging cable that has received a sharp blow, been dropped, run over, or damaged in any way.

8. **BATTERY CHEMICALS CAUSE SERIOUS BURNS.** Never allow the internal battery to come into contact with the skin, eyes, or mouth. If a damaged battery leaks chemicals, use rubber or neoprene gloves to dispose of it. If skin is exposed to battery fluids, wash with soap and water then rinse with vinegar. If eyes are exposed to battery chemicals, immediately rinse with water for 20 minutes and seek medical attention. Remove and dispose of contaminated

clothing.

9. **DO NOT OPERATE UNDER BLANKET AND PILLOW.** Excessive heating can occur and cause fire, electric shock, or injury to persons.

10. **STORING THE Foam Roller.** Store in a cool, dry place. Do not store the Foam Roller where temperatures may exceed 40°C/104°F such as in direct sunlight, in a vehicle, or metal building during the summer.

11. **DISPOSING OF BATTERIES.** The internal lithium-ion batteries for the Foam Roller are more environmentally friendly than some other types of batteries (e.g, nickel-cadmium). Always dispose of the Foam Roller according to government, regional, and local regulations. Contact a recycling agency in your area for recycling locations. Even discharged batteries contain some energy. Before disposing, use electrical tape to cover the Foam Roller's 12V connector to prevent the battery from shorting, which could cause a fire or explosion.

12. **DO NOT DISASSEMBLE.** Disassembly or incorrect reassembly may result in the risk of electric shock, fire, or exposure to battery chemicals. The warranty will be void if the Foam Roller is disassembled or if any parts have been removed.

13. **SERVICE.** If the Foam Roller is not working properly, has received a sharp blow, has been dropped, damaged, left outdoors, or dropped into water, do not use. Do not attempt to repair or disassemble the Foam Roller, which may result in an electric shock or fire.

14. **DO NOT USE WHILE BATHING OR IN THE SHOWER, TUB, OR SINK.** Do not place or store the Foam Roller where it can fall or be pulled into a tub or sink. Do not place in or drop into water or other liquid. Do not reach for an appliance that has fallen into water. Unplug immediately.

15. **THERMAL LIMITER.** The Foam Roller employs an automatic reset thermal limiter that shuts off the Foam Roller to prevent overheating and fire.

This device can be used by children aged from 8 years and above, and persons with reduced physical, sensory or mental capabilities, or lack of experience and knowledge if they have been

given supervision or instruction concerning use of the device in a safe way and understand the hazards involved. Children shall not play with the device. Cleaning and user maintenance shall not be used without supervision. This device contains batteries that are only replaceable by authorized personnel. This device contains batteries that are non-replaceable.

LIMITED WARRANTY

For full warranty information, please visit FlexirRecovery.com/warranty. To request a copy of the warranty by mail, you may send a request to the following address:

Flexir Recovery
Attn: Warranty Department
251 Little Falls Drive
Wilmington, DE 19807

Please note, this is not a return address or a retail location. No products or packages will be accepted at this location.

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Flexir, LLC. The Masonry, 151-156 Thomas St Suite 208, Dublin, Ireland